

Acquiring Luminescence Spectra using the ARC –NCL Spectral Measurement System + Boxcar Integrator





Disclaimer

- Safety –the first !!! This presentation is not manual. It is just brief set of rule to remind procedure for simple measurements. You should read manual first.
- Notwithstanding any language to the contrary, nothing contained herein constitutes, nor is intended to constitute, an offer, inducement, promise, or contract of any kind. The data contained herein are for informational purposes only and are not represented to be error free. Any links to non-UAB information are provided as a courtesy. They are not intended to constitute, nor do they constitute, an endorsement of the linked materials by the University of Alabama at Birmingham.




Luminescence

is a process in which an excited material emits light (electromagnetic radiation)

For example, the following types of luminescence caused by different excitation processes could be classified as :

Chemoluminescence, is the emission of light as the result of a chemical reaction

Light production in fireflies is due to a type of chemical reaction called bioluminescence



Photoluminescence (PL) is a process in which a material absorbs photons (electromagnetic radiation) and then re-radiates photons

Banknote photoluminescence after excitation by UV light from flash lamp



Electroluminescence (EL) is an process in which a material emits light in response to an electric current passed through it, or to a strong electric field.

Infrared electroluminescence of photodiode in the remote is stimulated by electrical current





Parameters of the Luminescence

- ❑ Luminescence **spectra** shows how intensity of the luminescence depends on wavelength
- ❑ Luminescence **lifetime** refers to the average time the molecule/ion stays in its excited state before emitting a photon (or how long luminescence could be observed)
- ❑ The luminescence **quantum yield** gives the efficiency of the luminescence. It is defined as the ratio of the number of photons emitted to the number of photons absorbed.


**Goal of the Lab: Measuring of the Luminescence spectra using
ARC –NCL Spectral System**

Experiment Background

1. First of all, we need excitation source to transfer our sample into excited (high energy) state. Here we will consider only optical excitation. It means that we will use optical radiation to excite the sample.
2. Second, we need an optical system to collect the luminescence and direct it to the detector. Also we need to select luminescence from the excitation radiation. For these purposes we can use optical filter or monochromator which can help select radiation only at required wavelength.
3. We need to detect optical signal (convert intensity of the optical radiation into electrical signals) and then convert signal into digital format




Experimental setup (Principal Scheme)




Requirements to optical excitation

A wavelength of the optical excitation pulses should be within absorption band of the studied samples !!!


Available commercial tunable solid-state lasers



OPO Ekspla




OPO Spectra Physics



Typical curve of the Alexandrite laser

Optical Setup

1) Avoid reflection of the excitation radiation to the measurement system!!!



2) Acceptance angle of Spectra Pro300i is $\Theta = 1/4$. Therefore lens diameter should be $D > L\Theta$ (for $L_2 = L_1 = 2F$ configuration $D > F/2$)

How to Choose a Detector ?


A **optical detector should** convert luminescence of the sample into a electrical signal.

Therefore the major requirement to the optical detector : to be sensitive at the wavelength of the to the luminescence photons




Operating Ranges for Judson Technologies detectors

Detectivity vs Wavelength for J10D Series InSb







Example of Detectivity for J15D Series HgCdTe




Attention, these detectors should cooled by liquid Nitrogen !!!

Operating Ranges for ARC detectors



Two type of the electronic registration system

depends on the ratio of the signal duration (T_s) and period between pulses (T_r); duty cycle



➤ Lock-in Amplifier for signals when $T_s \sim T_r$

➤ Boxcar Integrator for signals when $T_s \ll T_r$


Experiments with Boxcar Integrator



A gated integrator (also called a boxcar integrator or averager) integrates an analytical signal over a fixed time window. In pulsed experiments the integrator gate is synchronized with the analytical signal by a trigger. This method increases the signal-to-noise ratio by recording the voltage only when the signal is present, and ignoring time periods when there is no signal and only noise.

Spectra Pro 300i

Optical design: Czerny – Turner



Focal length: 300 mm

Aperture ratio: $f / 4$

Grating size: 68 x 68 mm

Grating mount: triple-grating turret

Grating #1 300gr/mm ($\lambda_{\text{blaze}}=$) ($\lambda_{\text{max}}=5.6$)

Grating #2 600 gr/mm ($\lambda_{\text{blaze}}=$) ($\lambda_{\text{max}}=2.8$)

Grating #3 150 gr/mm ($\lambda_{\text{blaze}}=$) ($\lambda_{\text{max}}=11.2$)


Standard slits: adjustable
from 10 μm to 3 mm wide;

Linear dispersion (nm/mm)@500 nm :

11- Grating #1 300gr/mm

5- Grating #2 600 gr/mm

21- Grating #3 150 gr/mm






Spectra Pro 300i Connections



To Power
Supply unit

To Control Unit
(NCL)



Experimental Setup and Wiring

Spectra Pro 300i

Thorlabs Si Photo Detector can be used. Detector should be set to "see" scattered light, **NOT** under direct beam.

sample


Detector

Pulsed Laser

Sample Luminescence signal

Averaged Signal

Ch1



Triggering signal

gate



Power switch

Description of Processes Above

- 1) Pulsed laser excites the sample
- 2) Si detector detects excitation pulse and triggers both Boxcar and Oscilloscope
- 3) Monochromator splits the PL signal from the sample and directs the desirable wavelength onto the detector
- 4) Triggered oscilloscope reads and displays the signal from the detector
- 5) Triggered Boxcar outputs a Gate Signal.
- 6) This Gate Signal needs to be adjusted and shifted onto the Luminescence signal using the controller knobs (see next slide for detailed description of Boxcar operation)
- 7) Boxcar receives the Luminescence signal coinciding with the gate and averages it increasing the amplitude and quality of the signal
- 8) Averaged signal is sent to PC through the NCL and Luminescence Spectrum is obtained

NOTE: Before proceeding with the steps described in the next slide bring the monochromator to the expected Luminescence Wavelength, ex: for Cr:ZnSe set monochromator to ~2600nm (Look at page 21 B)

Operation of Boxcar



Triggering - After wiring every component as shown above and firing the laser power on the Boxcar. Triggering is set using the RATE and RATE ADJUST knobs until green LED (TRIG) starts blinking. The Trigger Signal is used by Boxcar to initiate Gate Signal.

Delay – The **Gate Signal** needs to be **positioned onto the Luminescence Signal** (both of these signals should be displayed on the oscilloscope display, you are just trying to visually overlap the signal like shown on the Oscilloscope's display above). This is done using SCALE and MULTIPLIER knobs(in delay section). Read the time scale on the Oscilloscope this is roughly the scale you should set on the SCALE knob. Rotating MULTIPLIER knob you should be able to shift the Gate left and right on the screen.

Width - The width of the Gate Signal can be adjusted using the SCALE and MULTIPLIER knobs(in width section)

Signal - There are many sources of noise in this set up. These noises can be eliminated by off-setting the noise signal. To do this connect LAST SAMPLE (from OUTPUTS) to either ANALOG METER or DIGITAL METER port (to the left of the POWER Switch) the signal will be displayed on the analog or LED read-outs above. Cover the crystal with cardboard blocking Luminescence from crystal. Using the SENSITIVITY and INPUT OFFSET bring the signal to a small positive reading(in other words you just zeroed the background noise). Now remove the cardboard and check if the sensitivity is appropriate, if signal is too low rotate the SENSITIVITY knob to a lower voltage.

Averaging – Using the SAMPLES knob you can set how many signals are taken to output 1 averaged signal. The **integration time** in the “SpectraSense” software (next slide) needs to be set to integer multiples of $\left(\frac{1}{\text{laser pulse freq.}}\right) \cdot (\text{\#of samples})$ to prevent integration of signals corresponding to later wavelengths. For example, if laser has 10 pulses per second and SAMPLES is set to 3 Integration time can be 0.3, 0.6, 0.9...etc.

NOTE: Some of our detectors produce and negative signal, in this case there are two switches on the back of Boxcar that need to be switched to “inverted” position.





Operation Procedure


1. Switch on :

Boxcar Integrator
Detector Power Supply (see manual)

2. Switch on Spectra Pro 300i :


A-Switch on PC
B-Switch on Spectrometer power supply
C- Switch on Controller Unit (NCL)

Do not switch Laser on !!!



Operation Procedure

3. Run “SepctraSense” Program




Operation Procedure

4. Select Hardware Status Menu



Operation Procedure

4. Select grating (#1,#2, or #3) –A
5. To Select required wavelength : input wavelength, in nm and press “GOTO” -B




- 6 . Switch laser on according to the laser manual (see safety instruction !!!)

Operation Procedure

6. For Spectra Measuring:


6.A- Select Accusation mode



Operation Procedure

6. For Spectra Measuring:

6.B- Select Scan of the Monochromator #1




The screenshot shows the SpectraSense software window. The 'Scan Type' dropdown menu is open, displaying the following options: Scan Mono 1 (selected), Scan Mono 2, Time Based, Synchronized, Mono 1 vs 2, Mono 2 vs 1, and none selected. A yellow arrow points to the 'Scan Mono 1' option. The 'From' field is set to 0.000 nm, the 'To' field is set to 0.000 nm, the 'Increment' is 0.000 nm, and the 'Integration Time' is 100 mSec. The 'Acquisition Mode' is set to 'File', and the 'Storage' is set to 'averaged'. The 'Dark Subtraction' checkbox is unchecked. The 'Acquisition Routine' is set to 'None', and the 'Operator' field is empty. The 'Comments' field is empty. The 'Save data as' field is empty. The 'Routine' buttons are 'Save', 'Save As', and 'View'. The 'Settings' button is 'Current', and the 'Acquire' button is 'Acquire'. The 'Cancel' button is 'Cancel'. The bottom status bar shows 'Hardware Status', 'Survey Mode', 'Acquisition', 'Live Data', and 'Post Processing'.

Operation Procedure

6. For Spectra Measuring:

6.D- Select Acquisition Mode – Chanel #1



Operation Procedure

7. For Spectra Measuring:

7.E- Select Scan Parameters

The screenshot shows the SpectraSense software window with the following parameters and controls:


- Scan Type:** Scan Mono 1 (dropdown)
- MonoChromator 1:**
 - From:** 2440.000 nm (E1 - Initial Wavelength)
 - To:** 2600.000 nm (E2 - Final Wavelength)
 - Increment:** 2.000 nm (E3 - Increment)
 - Integration Time:** 100 mSec. (E4 - Integration Time at each wavelength)
 - Cycles:** 1 (E5 - Number of scanning cycles)
 - Delay:** 0 h 0 m 0 s
- Acquisition:**
 - Mode:** ch1 (dropdown)
 - File:** No File
 - Storage:** averaged (dropdown)
 - ☐ Dark Subtraction
- Acquisition Routine:** (dropdown) with a **Browse** button.
- Operator:** (dropdown)
- Comments:** (text area)
- Save data as:** (text field) with a **Browse** button.
- Routine:** (Save, Save As, View buttons)
- Settings:** (Current, Acquire, Cancel buttons)

A large green arrow points to the **Acquire** button at the bottom right of the window.

- E1-Initial Wavelength
- E2- Final Wavelength
- E3 -Increment
- E4- Integration Time at each wavelength
- E5- Number of scanning cycles

8. Press Acquire to start measurements

Results will be shown in the “Live Data” window




9. To save Spectra in the ASCII code (Text File)

9.1 Select File>Save

9.2 Select “Experiment Data”

9.3 Press “OK”

9.4 Select Directory and File Name





Operation Procedure

10. After work done, shutdown :

- ☐ Laser (according manual)
- ☐ Detector Power Supply Unit
- ☐ Boxcar Integrator
- ☐ NCL controller
- ☐ Monochromator power supply unite